Lympho 24+ Spin Medium / PBMC 24+ Spin Medium

Order No.
Size

60-00093-10

100 ml

60-00093-11

250 ml

60-00093-12

500 ml

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Separation scheme for PBMC enrichment with Lympho 24+ Spin Medium

At first carefully layer the sample material on top of the Lympho 24+ Spin Medium. Avoid a mixing of the two phases. Alternatively it's possible to use **pluriMate** tubes. The barrier prevents a mixing of sample material and density solution. Furthermore the tube enables pouring off the enriched cells after the first centrifugation step. After the density centrifugation aspirate the upper layer (Plasma and dilution buffer). Afterwards transfer the mononuclear cell layer to a new conical tube and wash the cells.


Fig. 1: Preparation of PBMC with Lympho 24+ Spin Medium. The sample material will be overlaid on top of density gradient medium.


Fig. 2.: Layers after density gradient centrifugation. The PBMC can be found on top of the Lympho Spin Medium. Erythrocytes, Granulocytes and dead cells will pass the medium and can be found at the bottom of the tube.

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Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.