

## PBMC 24+ Spin Medium (100 ml)

**Order No.:**60-00093-10



### Description

PBMC Spin 24+ Spin is a ready-to-use, sterile density medium for the isolation of PBMC from 8-24 hour old sample material (peripheral blood, buffy coat) through single step density gradient centrifugation. In contrast to PBMC Spin Medium, the use of PBMC24+ reduce the time dependent contamination with granulocytes and erythrocytes which occurs in PBMC isolation of blood older than 8 hours. With PBMC24+ Spin high purity and a yield up to 80% even 24 hours after sampling are still possible. The enriched cells can be used for a wide range of downstream applications such as magnetic activated cell sorting (MACS) or stem cell research. MACS is a trademark of Miltenyi Biotec GmbH.

### Additional Information

|                      |  |
|----------------------|--|
| Size                 | 100 ml   |
| Delivery Time (days) | 1-2  |
| Storage Condition    | Room temperature, protect from light                         |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |
| pH Range             | 6.8 - 7.4  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.