

Pestle (25 pcs.)



Order No.:05-00040-03

Description

Pestle for tissue disruption and preparing single-cell suspensions.

Recommended for use with pluriStrainer combined with Connector Ring.

Application



Pestle works well with pluriStrainer and Connector Ring:

1. Close Luer-Lock on Connector Ring to prevent flow-through
2. Add sample on top
3. Disrupt sample with pestle
4. Open Luer-Lock to allow flow-through

pluriSelect Life Science

Deutscher Platz 5a, 04103 Leipzig, DE

Phone: +49 341 333 858 0 | Fax: +49 341 333 858 99

support@pluriselect.com

<https://www.pluriselect.com/>

Additional Information

| | |
|-------------------|------------------|
| Size | 25 pcs. |
| Storage Condition | Room Temperature |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.