

## Connector Ring (50 pcs., non-sterile)

Order No.:41-50000-50



### Descripción

The Connector Ring is a sterile device to support active filtration, cell detachment, cell lysis or cell stimulation in combination with a pluriStrainer and a 50ml standard (conical) centrifuge tube. pluriStrainer, Connector Ring and a 50ml tube can be combined to a unit. This unit allows for a controlled flow through the sieve. When the Luer-lock of the Connector Ring is closed, sample or liquid remains on the sieve. An open Luer-lock starts the separation or the filtration process.

### Aplicación

#### The Connector Ring enables:

- Detachment of the targets from pluriBeads directly on the pluriStrainer
- Lyses of cells to extract cell specific RNA, DNA or proteins directly on the pluriStrainer
- Controlled flow through the cell strainer with a pump device or pipette
- Preparation of a single cell suspension from tissue and other materials
- Fast filtration of large volumes with efficient recovery rate.

### Características



**Flow control**

Allows to control the rate of flow by opening or closing the Luer-Lock, e.g. for physical dissociation of primary tissue (brain, spleen etc.).



**Low pressure**

If you add a syringe to the Connector Ring, it is possible to force low pressure to support the straining of rough sample material while pulling the piston.

**Additional Information**

Tamaño	50 pcs.
Tiempo de entrega	1-2
Variación	50 pcs, non-sterile
Condición de envío	La temperatura de la habitación
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.