

pluriSpin® Human NK Cell Depletion



| Order No. | Size |
|-------------|-------|
| 19-05002-31 | 2 ml |
| 19-05002-32 | 10 ml |

Description

pluriSpin® Human NK Cell Depletion is designed to isolate pure, viable and untouched nk cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels NK cells, while the isolated human nk cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified nk cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Informations complémentaires

| | |
|----------------------------|---|
| Espèces | Humain |
| Types de cellules | Cellules NK |
| Échantillon de matériel | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Pureté | >90% |
| Les conditions de stockage | Température ambiante |
| Déclaration réglementaire | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.