

## Re-Strainer (In-Line-Strainer) 70 µm - Sterile

Order No.:43-75070-40



### Description

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

### Candidature

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

### Informations complémentaires

|                    |                   |
|--------------------|-------------------|
| Délai de livraison | 1-2               |
| Taille des mailles | 70 µm             |
| Emballage          | 4 devices per bag |
| Taille             | 20 pcs.           |

|                            |  |
|----------------------------|--|
| Couleur                    | transparent  |
| Contrôle des flux          | oui  |
| Mesh Fixation              | Injected in housing for strong hold  |
| Tissu de la matière        | PET (polyéthylène téréphtalate)  |
| Matériau du boîtier        | LD-PE (polyéthylène basse densité)   |
| Stérilité                  | Stérile  |
| Stabilité                  | Good resistance to non-oxidizing acids and bases, fats, most organic solvents                  |
| Conditions d'expédition    | Room Temperature   |
| Les conditions de stockage | Température ambiante   |
| Déclaration réglementaire  | For research use only  |
| Informations juridiques    | <b>The product is for research and development only, not for diagnostic or therapeutic use</b> |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.