

PBMC 24+ Spin Medium (250 ml)

Order No.:60-00093-11



Description

PBMC Spin 24+ Spin is a ready-to-use, sterile density medium for the isolation of PBMC from 8-24 hour old sample material (peripheral blood, boffy coat) through single step density gradient centrifugation. In contrast to PBMC Spin Medium, the use of PBMC24+ reduce the time dependend contamination with granulocytes and erythrocytes which occurs in PBMC isolation of blood older than 8 hours. With PBMC24+ Spin high purity and a yield up to 80% even 24 hours after sampling are still possible. The enriched cells can be used for a wide range of downstream applications such as magnetic activated cell sorting (MACS) or stem cell research. MACS is a trademark of Miltenyi Biotec GmbH.

Additional Information

Size	250 ml
Delivery Time (days)	1-2
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 - 7.4

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.