

# Steel Basket-Strainer 500 µm (10 pcs., non-sterile)

Order No.:43-76500-50



#### **Descrizione**

 $500 \ \mu m$  Steel Basket Strainers are stainless steel mesh filters that consist of individual strands woven into a mesh. The Strainers are characterized by precise mesh openings, percent open area and mesh thickness.

Steel Basket-Strainer are available with five different pore sizes 20  $\mu$ m, 50  $\mu$ m, 100  $\mu$ m, 200  $\mu$ m and 500  $\mu$ m.

Steel Basket Strainer consits of a SS 304 stainless steel ring and a SS 316 stainless steel mesh, for the best efficiency of your application. The steel mesh filter can be reused. Therefore the strainers can be regenerated by ultrasonic bath or backwashing the system. After regeneration, the steel filter will be ready for use again. Additionally heat sterilization or autoclaving is possible. The stainless steel filters are very robust and can tolerate high temperatures and high pressure. These steel filters are usable for liquids or compressed air. Due to the strenth of the mesh offer an absolute retention rating in a wide range of filtration applications.

#### **Applicazione**

- Extract liquids from porous sample material
- · Dissociation of tissue and cells
- Preparation of single cell suspension
- Alternative strainer to gauze filtration
- Alternative strainer to plastic cell strainer
- Usable at higher temperatures
- Usable for petroleum and chemical industries, in the aerospace industry or in food processing
- Particle retention in cosmetic
- Water treatment
- Syrup filtration

### **Additional Information**

Tempi di consegna	1-2
Dimensione della maglia	500 μm
Variazione	10 pcs, Stainless Steel Mesh
Dimensione	10 pcs
Materiale del tessuto	Steel
Mesh Type	woven
Sterilità	Non-Sterile
Centrifugabile	2000 x g (max.)
Dissociazione dei tessuti	sì
Disinfettabile	sì
Condizioni di spedizione	Temperatura ambiente
Condizione di conservazione	Temperatura ambiente
Dichiarazione di regolamentazione	For research use only. Not for use in diagnostic procedures.

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or theurapeutic use.