

pluriSpin® Human NK Cell Depletion

411	Order No.	Size Size
15	19-05002-31	2 ml
pluriSpin nativo Celi Becantilla	19-05002-32	10 ml
Tenan		

概要

pluriSpin® Human NK Cell Depletion is designed to isolate pure, viable and untouched nk cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels NK cells, while the isolated human nk cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified nk cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

しゅるい	人間
細胞型	NK細胞
サンプル材料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
純潔	>90%
保管状態	室温
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.