

Detachment Buffer (Buffer C & Buffer D)

Order No.:60-00046-12



概要

Detachment Buffer - includes Buffer C (Detachment Activation Buffer, 60-00045-12) and Buffer D (Detachment Concentrate, 60-00040-12). The product is used for the detaching of the pluriBeads after cell separation.

Additional Information

大きさ	10 * 1.8 ml + 10 ml
保管状態	Buffer C 4-8 °C; Buffer D -20°C
きせいぶん	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.