

## Pestle (25 pcs.)



**Order No.:**05-00040-03

### Description

Pestle for tissue disruption and preparing single-cell suspensions.

Recommended for use with pluriStrainer combined with Connector Ring.

### Application



Pestle works well with pluriStrainer and Connector Ring:

1. Close Luer-Lock on Connector Ring to prevent flow-through
2. Add sample on top
3. Disrupt sample with pestle
4. Open Luer-Lock to allow flow-through

**pluriSelect Life Science**

Deutscher Platz 5a, 04103 Leipzig, DE

Phone: +49 341 333 858 0 | Fax: +49 341 333 858 99

support@pluriselect.com

<https://www.pluriselect.com/>

### **Additional Information**

Size	25 pcs.
Storage Condition	Room Temperature

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.