

PBMC 24+ Spin Medium (250 ml)

Order No.:60-00093-11



Description

PBMC Spin 24+ Spin is a ready-to-use, sterile density medium for the isolation of PBMC from 8-24 hour old sample material (peripheral blood, buffy coat) through single step density gradient centrifugation. In contrast to PBMC Spin Medium, the use of PBMC24+ reduce the time dependent contamination with granulocytes and erythrocytes which occurs in PBMC isolation of blood older than 8 hours. With PBMC24+ Spin high purity and a yield up to 80% even 24 hours after sampling are still possible. The enriched cells can be used for a wide range of downstream applications such as magnetic activated cell sorting (MACS) or stem cell research. MACS is a trademark of Miltenyi Biotec GmbH.

Additional Information

| | |
|----------------------|--|
| Size | 250 ml |
| Delivery Time (days) | 1-2 |
| Storage Condition | Room temperature, protect from light |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |
| pH Range | 6.8 - 7.4 |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.