

## PBMC 24+ Spin Medium (100 ml)

**Order No.:** 60-00093-10



### Description

PBMC Spin 24+ Spin is a ready-to-use, sterile density medium for the isolation of PBMC from 8-24 hour old sample material (peripheral blood, buffy coat) through single step density gradient centrifugation. In contrast to PBMC Spin Medium, the use of PBMC24+ reduce the time dependent contamination with granulocytes and erythrocytes which occurs in PBMC isolation of blood older than 8 hours. With PBMC24+ Spin high purity and a yield up to 80% even 24 hours after sampling are still possible. The enriched cells can be used for a wide range of downstream applications such as magnetic activated cell sorting (MACS) or stem cell research. MACS is a trademark of Miltenyi Biotec GmbH.

### Additional Information

Size	100 ml
Delivery Time (days)	1-2
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 – 7.4

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.