

## pluriSpin® Human NK Cell Depletion

	Order No.	Size
	19-05002-31	2 ml
pluriSpin nanve con Separation	19-05002-32	10 ml
A Manual A		

## **Description**

pluriSpin® Human NK Cell Depletion is designed to isolate pure, viable and untouched nk cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels NK cells, while the isolated human nk cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified nk cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

## **Additional Information**

Species	Human	
Cell Types	NK cells	
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution	
Purity	>90%	
Storage Condition	Room Temperature	
Regulatory Statement	For research use only. Not for use in diagnostic procedures,	

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or theurapeutic use.