

Pestle (25 pcs.)



Order No.:05-00040-03

Description

Pestle for tissue disruption and preparing single-cell suspensions.

Recommended for use with pluriStrainer combined with Connector Ring.

Application



Pestle works well with pluriStrainer and Connector Ring:

1. Close Luer-Lock on Connector Ring to prevent flow-through
2. Add sample on top
3. Disrupt sample with pestle
4. Open Luer-Lock to allow flow-through

PluriSelect USA

2030 Gillespie Way, Unit 103, CA 92020 El Cajon, US
Phone: +1 619-202-4297 | Fax: +1 619-828-8003
sales-support@pluriselect-usa.com
<https://www.pluriselect.com/>

Additional Information

Size	25 pcs.
Storage Condition	Room Temperature

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.