

EDTA Solution (0.2M in PBS)



| Order No. | Size | Sterility |
|-------------|--------|-----------|
| 60-00030-10 | 10 ml | Sterile |
| 60-00030-11 | 50 ml | Sterile |
| 60-00030-12 | 100 ml | Sterile |

Description

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as a chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

Additional Information

| | |
|----------------------|--|
| Area of interest | Cell Biology |
| Sterility | Sterile |
| pH Range | pH 7.4 |
| Shipping Condition | Room Temperature |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.