

## universal M-pluriBead® anti-hu

**Order No.:**31-RBaHU-20



### Description

pluriBead® universal are developed for the use with external antibodies. They are coated with a polyclonal antibody produced from the serum of a host animal. (2ml)

### Additional Information

|                      |                                                                               |
|----------------------|-------------------------------------------------------------------------------|
| Size                 | 2 ml                                                                          |
| Species              | Human                                                                         |
| Cell Types           | Any                                                                           |
| Sample Material      | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity               | >90%                                                                          |
| Storage Condition    | 4-8 °C                                                                        |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.