

CD43 M-pluriBead® anti-hu

Order No.:19-04300-20



Description

CD43 pluriBead® suspension for the positive cell isolation of CD43+ target cells from any human sample material. (2ml)

Additional Information

Size	2 ml
Species	Human
Cell Types	T cells, Granulocytes
Sample Material	Whole Blood, PBMC, Tissue, Buffy Coat, Cord Blood, Bone Marrow
Purity	>90%
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.