

## S-pluriBead® Maxi Reagent Kit

**Order No.:**70-50010-12



### Description

Maxi S-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer, detachment activation buffer, 10x 1,8ml detachment buffer concentrate, 10 S-pluriStrainer, 10 Funnel, 10 connector rings

### Additional Information

Size	10
Delivery Time (days)	1-2
Species	Human
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	Room Temperature; Buffer 4-8 °C; Buffer D -20°C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.