

10% BSA Stock Solution in 1x PBS



| Order No. | Size | Sterility |
|-------------|--------|-----------|
| 60-00020-10 | 50 ml | Sterile |
| 60-00020-11 | 250 ml | Sterile |
| 60-00020-12 | 500 ml | Sterile |

Description

10% BSA stock solution in 1x PBS can used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

Additional Information

| | |
|----------------------|--|
| Delivery Time (days) | 1-2 |
| Area of interest | Cell Biology |
| Sterility | Sterile |
| pH Range | pH 7.4 |
| Shipping Condition | Room Temperature |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

