

Syringe Strainer 5 µm (12 pcs., non-sterile)

Order No.:43-71005-50



Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

| | |
|----------------------|--|
| Delivery Time (days) | 1-2 |
| Mesh Size | 5 µm |
| Variation | 12 pcs, non-sterile |
| Size | 12 pcs. |
| Flow Control | Yes |
| Housing Material | LD-PE (Low Density Polyethylen) |
| Mesh Type | woven |
| Sterility | Non-Sterile |
| Shipping Condition | Room Temperature |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.