

## Syringe Strainer 40 µm (12 pcs., non-sterile)



**Order No.:**43-71040-50

### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

Mesh Size	40 µm
Variation	12 pcs, non-sterile
Size	12 pcs.
Flow Control	Yes
Housing Material	LD-PE (Low Density Polyethylen)
Mesh Type	woven
Sterility	Non-Sterile
Shipping Condition	Room Temperature
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.