

PLT Spin Medium (500 ml)

Order No.: 60-00094-12



Description

PLT Spin Medium is a ready-to-use, sterile density medium for the isolation of platelets in high yield from fresh human peripheral blood through single step centrifugation.

Additional Information

| | |
|----------------------|--|
| Size | 500 ml |
| Delivery Time (days) | 1-2 |
| Storage Condition | Room temperature, protect from light |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |
| pH Range | 6.8 – 7.4 |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.