

## pluriSpin® Human Monocyte Enrichment



| Order No.   | Size  |
|-------------|-------|
| 19-01001-31 | 2 ml  |
| 19-01001-32 | 10 ml |

### Description

pluriSpin® Human Monocyte Enrichment is designed to isolate pure, viable and untouched monocytes from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, B cells, NK cells and granulocytes, while the negatively isolated human monocytes are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

|                      |   |
|----------------------|---|
| Species              | Human   |
| Cell Types           | Monocytes   |
| Sample Material      | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity               | >90%  |
| Storage Condition    | Room Temperature  |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.