

1x PBS Solution (10 * 500 ml)

Order No.:60-00000-12



描述

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

Components

- NaCl (136,89 mMol/L)
- Na₂ HPO₄ • 2H₂ O (10mMol/L)
- KH₂ PO₄ (1,98 mMol/L)
- KCl (2,68 mMol/L)

Additional Information

大小	10 * 500 ml
交付时间	1-2
感兴趣的领域	细胞生物学
不育症	无菌型
相位范围	pH 7.4
运输条件	室温

储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.