

CD193 M-pluriBead® anti-hu

Order No.:19-19300-20



描述

CD193 pluriBead® suspension for the positive cell isolation of CD193+ target cells from any human sample material. (2ml)

Additional Information

大小	2 ml
品种	人类
细胞类型	嗜碱性粒细胞, 嗜酸性粒细胞 粒细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.