

## CD16 M-pluriBead® anti-hu



**Order No.:**19-01600-20

### 描述

CD16 pluriBead® suspension for the positive cell isolation of CD16+ target cells from any human sample material. (2ml)

### Additional Information

|       |   |
|-------|---|
| 大小    | 2 ml  |
| 品种    | 人类  |
| 细胞类型  | 嗜中性粒细胞  |
| 样品料   | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 纯度    | >90%  |
| 储存条件  | 4-8 °C  |
| 规范性声明 | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.