

## CD11c S-pluriBead® anti-ms

**Order No.:**29-011c0-10



### 描述

CD11c pluriBead® suspension for the positive cell isolation of CD11c+ target cells from any murine sample material. (2ml)

### Additional Information

大小	2 ml
品种	小鼠(鼠类)
细胞类型	活化的T细胞, 树状细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.