

F4/80 S-pluriBead® anti-ms

Order No.:29-F4800-10



描述

F4/80 pluriBead® suspension for the positive cell isolation of F4/80+ target cells from any murine sample material. (2ml)

Additional Information

| | |
|-------|---|
| 大小 | 2 ml |
| 品种 | 小鼠(鼠类) |
| 细胞类型 | 单核细胞, 巨噬细胞 |
| 样品料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 纯度 | >90% |
| 储存条件 | 4-8 °C |
| 规范性声明 | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.