

## CD271 S-pluriBead® anti-hu

**Order No.:**19-27100-10



### 描述

CD271 pluriBead® suspension for the positive cell isolation of CD271+ target cells from any human sample material. (2ml)

### Additional Information

大小	2 ml
品种	人类
细胞类型	间质细胞, 许旺细胞, 干/基质细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.