

10% BSA Stock Solution in 1x PBS



Order No.	Size	Sterility
60-00020-10	50 ml	无菌型
60-00020-11	250 ml	无菌型
60-00020-12	500 ml	无菌型

描述

10% BSA stock solution in 1x PBS can used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

Additional Information

交付时间	1-2
感兴趣的领域	细胞生物学
不育症	无菌型
相位范围	pH 7.4
运输条件	室温
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

