

EDTA Solution (0.2M in PBS)



Order No.	Size	Sterility
60-00030-10	10 ml	无菌型
60-00030-11	50 ml	无菌型
60-00030-12	100 ml	无菌型

描述

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as a chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

Additional Information

感兴趣的领域	细胞生物学
不育症	无菌型
相位范围	pH 7.4
运输条件	室温
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.