

## pluriSpin® Human CD4<sup>+</sup> T Cell Enrichment


**Order No.**
**Size**

19-03001-31

2 ml

19-03001-32

10 ml

**描述**

pluriSpin® Human CD4<sup>+</sup> T Cell Enrichment is designed to isolate pure, viable and untouched CD4<sup>+</sup> T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, cytotoxic T cells, B cells, NK cells and granulocytes, while the negatively isolated human CD4<sup>+</sup> T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD4<sup>+</sup> T cells are present as a highly enriched population at the interface between the plasma and the density medium.

**Additional Information**

品种	人类
细胞类型	辅助T细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures,

**Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.