

pluriSpin® Human CD8⁺ Cell Depletion


Order No.
Size

19-04002-31

2 ml

19-04002-32

10 ml

描述

pluriSpin® Human CD8⁺ Cell Depletion is designed to isolate pure, viable and untouched cytotoxic t cell-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Cytotoxic T cells, while the isolated human cytotoxic t cell-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified cytotoxic t cell-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|-------|---|
| 品种 | 人类 |
| 细胞类型 | 细胞毒性T细胞 |
| 样品料 | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 纯度 | >90% |
| 储存条件 | 室温 |
| 规范性声明 | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.