

pluriSpin® Human NK Cell Enrichment



Order No.	Size
19-05001-31	2 ml
19-05001-32	10 ml

描述

pluriSpin® Human NK Cell Enrichment is designed to isolate pure, viable and untouched NK cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, T cells, B cells and granulocytes, while the negatively isolated human NK cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified NK cells are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

品种	人类
细胞类型	NK细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>80%
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.