

pluriSpin® Human PBMC Depletion (10 ml)

Order No.:19-08002-32



描述

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

Additional Information

大小	10 ml
品种	人类
细胞类型	非PBMC细胞
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.