

## ÜberStrainer 1 µm



**Order No.:**43-70001-03

### 描述

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

### Additional Information

网眼尺寸	1 µm
大小	10 pcs
布料	PET(聚乙二醇酯)
外壳材料	PP (聚丙烯)
Mesh Type	woven
不育症	无菌型
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

