

## Syringe Strainer 300 µm (12 pcs., sterile)

**Order No.:**43-71300-51



### 描述

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

网眼尺寸	300 µm
变化	12 pcs sterile
大小	12 pcs.
流量控制	是
外壳材料	LD-PE (低密度聚乙烯)
Mesh Type	woven
不育症	无菌型
运输条件	室温
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.