

Syringe Strainer 15 µm



Order No.	Sterility	Variant
43-71015-50	非无菌型	12 pcs, non-sterile
43-71015-51	无菌型	12 pcs sterile

描述

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

网眼尺寸	15 µm
大小	12 pcs.
流量控制	是
外壳材料	LD-PE (低密度聚乙烯)
Mesh Type	woven
不育症	无菌型
运输条件	室温
储存条件	室温
规范性声明	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

