

Re-Strainer (In-Line-Strainer) 1 µm - Sterile

Order No.:43-75001-40



描述

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

应用

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

Additional Information

交付时间	1-2
网眼尺寸	1 µm
变化	20 pcs, sterile
包装	4 devices per bag

大小	20 pcs.
颜色	红色
流量控制	是
Mesh Fixation	Injected in housing for strong hold
布料	PET(聚乙二醇酯)
外壳材料	LD-PE (低密度聚乙烯)
不育症	无菌型
稳定性	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
运输条件	Room Temperature
储存条件	室温
规范性声明	For research use only
法律信息	The product is for research and development only, not for diagnostic or therapeutic use

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.