

## Re-Strainer (In-Line-Strainer) 100 µm - Sterile

**Order No.:**43-75100-40



### 描述

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

### 应用

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

### Additional Information

|      |                   |
|------|-------------------|
| 交付时间 | 1-2               |
| 网眼尺寸 | 100 µm            |
| 包装   | 4 devices per bag |
| 大小   | 20 pcs.           |

|               |  |
|---------------|--|
| 颜色            | 黄  |
| 流量控制          | 是  |
| Mesh Fixation | Injected in housing for strong hold  |
| 布料            | PET(聚乙二醇酯)   |
| 外壳材料          | LD-PE (低密度聚乙烯)   |
| 不育症           | 无菌型  |
| 稳定性           | Good resistance to non-oxidizing acids and bases, fats, most organic solvents                  |
| 运输条件          | Room Temperature   |
| 储存条件          | 室温   |
| 规范性声明         | For research use only  |
| 法律信息          | <b>The product is for research and development only, not for diagnostic or therapeutic use</b> |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.