

## pluriSpin® Human CD8<sup>+</sup> T Cell Enrichment (2 ml)

**Order No.:**19-04001-31



### 描述

pluriSpin® Human CD8<sup>+</sup> T Cell Enrichment is designed to isolate pure, viable and untouched cytotoxic T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, helper T cells, B cells, NK cells and granulocytes, while the negatively isolated human cytotoxic T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified cytotoxic T cells are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

大小	2 ml
品种	人类
细胞类型	细胞毒性T细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.

