

pluriSpin® Human Monocyte Enrichment (2 ml)

Order No.:19-01001-31



描述

pluriSpin® Human Monocyte Enrichment is designed to isolate pure, viable and untouched monocytes from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, B cells, NK cells and granulocytes, while the negatively isolated human monocytes are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

大小	2 ml
品种	人类
细胞类型	单核细胞
样品料	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
纯度	>90%
储存条件	4-8 °C
规范性声明	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

