

## pluriSpin® Human CD19<sup>+</sup> Cell Depletion (2 ml)

**Order No.:**19-06002-31



### 描述

pluriSpin® Human CD19<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched b cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels B cells, while the isolated human b cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified b cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

|       |   |
|-------|---|
| 大小    | 2 ml  |
| 品种    | 人类  |
| 细胞类型  | B细胞   |
| 样品料   | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| 纯度    | >90%  |
| 储存条件  | 4-8 °C  |
| 规范性声明 | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

