

High Density Spin Medium 500 ml (1.1 g/ml)

Order No.:60-01103-12



Description

Density Diluent Medium is a ready to use, sterile medium for the preparation of different density working solutions for the enrichment of cells, bacteria and other particles. The working solutions can be prepared by the dilution of High Density Spin Medium with Density Diluent Medium.

Additional Information

| | |
|----------------------|--|
| Size | 500 ml |
| Delivery Time (days) | 1-2 |
| Shipping Condition | Room Temperature |
| Storage Condition | Room temperature, protect from light |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |
| pH Range | 6.8 – 7.4 |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.