

CD14 S-pluriBead® anti-ms

Order No.:29-01400-10



Description

CD14 pluriBead® suspension for the positive cell isolation of CD14+ target cells from any murine sample material. (2ml)

Additional Information

| | |
|----------------------|-------------------------------------------------------------------------------|
| Size | 2 ml |
| Species | Mouse (Murine) |
| Cell Types | Monocytes |
| Sample Material | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity | >90% |
| Storage Condition | 4-8 °C |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.