

## pluriSpin® Human CD45<sup>+</sup> Cell Depletion





Order No.	Size
19-09002-31	2 ml
19-09002-32	10 ml

### Description

pluriSpin® Human CD45<sup>+</sup> Cell Depletion is designed to isolate pure, viable and untouched CD45<sup>+</sup>-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45<sup>+</sup> cells, while the isolated population free of CD45<sup>+</sup> are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45<sup>+</sup> cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Species	Human
Cell Types	CD45- Cell
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.