

10x PBS Solution (500 ml)

Order No.:60-00010-11



Beschreibung

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- RNase/DNase free
- Protease free

Components

- NaCl (136,89 mMol/L)
- Na₂HPO₄ • 2H₂O (10mMol/L)
- KH₂PO₄ (1,98 mMol/L)
- KCl (2,68 mMol/L)

Additional Information

Lieferumfang	500 ml
Sterilität	steril
ph-Bereich	pH 6.8
Lieferbedingungen	Raumtemperatur
Lagerbedingungen	Raumtemperatur

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.