

## 1x PBS Solution (500 ml)



**Order No.:**60-00000-10

### Beschreibung

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

### Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub>HPO<sub>4</sub> • 2H<sub>2</sub>O (10mMol/L)
- KH<sub>2</sub>PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Additional Information

Lieferumfang	500 ml
Lieferzeit	1-2
Interessengebiet	Zellbiologie
Sterilität	steril
ph-Bereich	pH 7.4
Lieferbedingungen	Raumtemperatur

Lagerbedingungen	Raumtemperatur
Regulierungserklärung	For research use only. Not for use in diagnostic procedures,

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.