


## 10% BSA stock solution in 1x PBS (500 ml)

**Order No.:**60-00020-12



### Beschreibung

10% BSA stock solution in 1x PBS can be used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

### Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

### Additional Information

Lieferumfang	500 ml
Lieferzeit	1-2
Interessengebiet	Zellbiologie
Sterilität	steril
ph-Bereich	pH 7.4
Lieferbedingungen	Raumtemperatur
Lagerbedingungen	4-8 °C
Regulierungserklärung	For research use only. Not for use in diagnostic procedures,

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.